

# caché™

nipple shields


Nipple hypersensitivity and projection are the often unavoidable consequence of breast enhancement procedures. As a result, patients frequently experience weeks or months of protracted post-operative nipple discomfort. Abrasion from surgical brassieres or wraps aggravates the nipple sensitivity, escalating the discomfort into persistent pain. Prominent nipple projection and fluid discharge commonly accompany the irritation, leading to patient embarrassment and the reluctant concealment of their new surgical results.

Caché™ shields have been designed to simply and inexpensively resolve the dilemma.

Caché™ shields are reusable, self-adhering nipple shields and covers. Constructed of soft, hypoallergenic silicone, Caché™ shields are flexible, comfortable and inconspicuous.

The interior surface features a patent-pending nipple pocket that houses and protects the nipple from any material abrasion. In addition to protection, this pocket masks nipple visibility through even the most thin or sheer garments.

Once applied and warmed to body temperature, translucent Caché™ shields become fully transparent. This functional invisibility allows Caché™ to adapt to any skin coloring.



© Copyright 2008 Becon Medical Limited 44BM0119Rev.A

becon

# caché™

nipple shields



A single set of Caché™ shields can be rinsed and reused for countless applications. Once hypersensitivity and projection diminish, the patient can continue to use their Caché™ shields for nipple concealment.

Physicians typically counsel their prospective patients about the benefits of Caché™ shields during pre-operative counseling. They may also furnish them during post-operative examination to satisfy the patient's complaint of hypersensitivity and projection.

Caché™ offers a simple and inexpensive solution to a common and needless problem.

*Patents pending*

Becon Medical Limited  
3530 East Grant Road  
Tucson, AZ 85716  
520-881-0480  
866-431-0480  
Fax: 520-881-0941  
[www.beconmedical.com](http://www.beconmedical.com)

<b>Catalog Number</b>	<b>Case Quantity</b>
BE-500	Box of two
BE-501	Case of four boxes (two per box)
BE-502	Case of ten boxes (two per box)

© Copyright 2008 Becon Medical Limited  
44BM0119Rev.A