


# Brecise™

precision measuring system for surgical breast procedures



## Brecise™... when you need to be precise

As an aesthetic surgeon, your goal of a perfect result for every breast procedure depends on exact pre and intra-operative measurement to guide patient skin marking. The rulers, tapes, strings and sutures currently employed to identify those measurements are unable to pinpoint these critical locations relative to fixed anatomical features.

Brecise™, a new and unique breast measuring device, has been created to easily and accurately quantify and document exact topographical locations critical to the conduct of consistent breast surgical procedures. Measurements are given in both distance and degree of angle, for best results in *all* major breast surgeries.

## Comprehensive

- Ensures accurate nipple-areolar complex positioning for breast reduction, breast reconstruction, mastopexy or any other procedure where nipple removal or movement is required.
- Identifies and documents pre-operative nipple positions and asymmetries for breast augmentations.
- \* Expedites pre and post operative patient skin marking.

## Inexpensive

- Brecise™ is packaged and economically priced as a single use, sterile, disposable product.


© Copyright Becon Medical Limited 44BM007Rev.A

becon


# Brecise™

precision measuring system for surgical breast procedures


**Sternal Notch to Nipple**




**Sternal Notch to Clavicle**




**Clavicle to Nipple**



**Inframammary Crease Angle**



**Nipple to Nipple**



## Versatile

Brecise™ enables multiple measurement options for comprehensive breast procedure skin marking and documentation:

- Sternal notch to nipple
- Sternal notch to clavicle
- Clavicle to nipple
- Inframammary crease angle
- Nipple to nipple

## Convenient

Included in each package is a visual data entry log for all Brecise™ measurements in both distance and angles, ready for insertion in the patient's chart.

Catalog Number	Case Quantity
BE-100	10 units per case
BE-100R	5 units per case

*Multiple U.S. and International patents pending*

Becon Medical Limited  
4640 East Sunrise Drive  
Tucson, AZ 85718  
520-881-0480  
866-431-0480  
Fax: 520-881-0941  
[www.beconmedical.com](http://www.beconmedical.com)

© Copyright Becon Medical Limited  
44BM007Rev.A