

Building your practice through patient safety



LipoStructure®

disposable instrumentation for microinfiltration

LipoStructure®* is the technique for microinjection procedures pioneered and popularized by Sydney Coleman, M.D. Becon Medical Limited has introduced a complete portfolio of advanced, disposable LipoStructure® instruments developed for these exacting procedures, which enable the surgeon to gently harvest and transplant intact cells with a minimum of trauma and a high probability of tissue survival.

DISPOSABLE CANNULAE - THE ONLY CHOICE FOR PROCEDURAL INTEGRITY AND SAFETY

After repeated autoclaving, reusable cannulae will corrode and pit, both internally and externally. Blood and tissue residues harden and accumulate over time, clogging the interior lumen. The combined effect of pitting and internal debris negates the intent of infiltration cannulae, which are designed to preserve the viability of fragile cells. Adopting single use LipoStructure® cannulae assures the smooth transfer of aspirant and enhances the probability of a successful graft.

Recent concern about cross-contamination has focused international attention on the obstacles to successfully sterilizing surgical cannulae. Cannulae, with their closed tips and thin bore, present obvious impediments to thorough cleaning. Since an instrument must be completely clean to sterilize, the inevitable presence of residual tissue, blood, fluids and other contaminants in the lumen of a cannula casts doubt on the efficacy of their sterilization.


LipoStructure® disposable cannulae allay all concerns of contamination from prior use, and will distinguish your practice as a unique provider of safe and effective microinfiltration instruments.

* LipoStructure® is a registered trademark of Sydney Coleman, M.D.


- **Sterile, disposable instrumentation eliminates danger of cross-contamination.**
- **Color coded for easy identification.**
- **Ultra thin wall for maximum internal flow.**
- **NEW Hood tip style included in this series.**
- **NEW micro cannulae.**
- **NEW facial and body dissectors.**




NEW LipoStructure® Body Cannula set, 15cm and 23cm




NEW LipoStructure® Mini Cannula set, Styles I and III with obturator



LipoStructure® Cannulae, Tip Styles I, II and III, 7cm and 9cm




Becon Infiltration Cannulae, 15cm - 35cm



NEW LipoStructure® W Dissector, 9cm Facial and 15cm Body



NEW LipoStructure® Hood Cannulae, 15cm and 23cm



LipoStructure® Aspiration Cannulae, 15cm and 23cm



LipoStructure® Style I Tip
LipoStructure® Style II Tip
LipoStructure® Style III Tip

Becon Medical Limited
3530 East Grant Road
Tucson, AZ 85716
520-881-0480
866-431-0480
Fax: 520-881-0941

© Copyright Becon Medical Limited 44BM0113Rev.A

www.beconmedical.com

Cat. No.	Description	Color
LS-DI7	LipoStructure Cannula, Style I, 7cm	Red
LS-DI9	LipoStructure Cannula, Style I, 9cm	Red
LS-DII7	LipoStructure Cannula, Style II, 7cm	White
LS-DII9	LipoStructure Cannula, Style II, 9cm	White
LS-DII9C	LipoStructure Cannula, Style II, Concave, 9cm	White
LS-DIII7	LipoStructure Cannula, Style III, 7cm	Blue
LS-DIII9	LipoStructure Cannula, Style III, 9cm	Blue
LS-DIII9C	LipoStructure Cannula, Style III, Concave, 9cm	Blue
LS-DI3	LipoStructure Cannula, Mini, Style I, 3cm with Obturator	Red Green
LS-DI6	LipoStructure Cannula, Mini, Style I, 6cm with Obturator	Red Yellow
LS-DIII3	LipoStructure Cannula, Mini, Style III, 3cm with Obturator	Blue Green
LS-DIII6	LipoStructure Cannula, Mini, Style III, 6cm with Obturator	Blue Yellow
LS-DW9	LipoStructure Facial W Dissector, 9cm	White
LS-DW15	LipoStructure Body W Dissector, 15cm	White
LS-DASP15	LipoStructure Aspiration Cannula, 15cm	Black
LS-DASP23	LipoStructure Aspiration Cannula, 23cm	Black
LS-DASP23C	LipoStructure Aspiration Cannula, Curved 23cm	Black

Cat. No.	Description	Color
LS-DH1015	LipoStructure Hood Cannula, 10ga x 15cm	Orange
LS-DH1023	LipoStructure Hood Cannula, 10ga x 23cm	Orange
LS-DH1215	LipoStructure Hood Cannula, 12ga x 15cm	Green
LS-DH1223	LipoStructure Hood Cannula, 12ga x 23cm	Green
LS-DB1215	LipoStructure Body Cannula, Style I, 12ga x 15cm	Red
LS-DB1223	LipoStructure Body Cannula, Style I, 12ga x 23cm	Red
LS-DB1415	LipoStructure Body Cannula, Style I, 14ga x 15cm	Red
LS-DB1423	LipoStructure Body Cannula, Style I, 14ga x 23cm	Red
OTHER INFILTRATION PRODUCTS FROM BECON:		
BE-1215	Infiltration Cannula, 12ga x 15cm	Blue
BE-1223	Infiltration Cannula, 12ga x 23cm	Blue
BE-1230	Infiltration Cannula, 12ga x 30cm	Blue
BE-1235	Infiltration Cannula, 12ga x 35cm	Blue
BE-1415	Infiltration Cannula, 14ga x 15cm	White
BE-1423	Infiltration Cannula, 14ga x 23cm	White
BE-1430	Infiltration Cannula, 14ga x 30cm	White
BE-1615	Infiltration Cannula, 16ga x 15cm	Yellow
BE-1620	Infiltration Cannula, 16ga x 20cm	Yellow
BE-LL3	Luer to Luer Connector	

**INFILTRATION CANNULAE ARE PACKAGED 10/BOX
LUER TO LUER CONNECTORS ARE PACKAGED 20/BOX
ALL OTHER PRODUCTS ARE PACKAGED 5/BOX**